



Stallion Medical Inc.

Face Mask

3-Layer Disposable
(box w/50pcs)



Specification

Material: Non woven and meltblown grade 95

Size: 175 x 95mm

Layers: 3 Layers

FE (filter efficiency %): $\geq 90\%$


Airflow Resistance [PA]: Ventilation resistance for gas exchange on both sides ≤ 15

Executive Standard: GB/T32610-2016

Storage: Temperature: $(25\pm 5)^\circ\text{C}$ Humidity: $(30\pm 10)\%$

Validity: 2 Years

3 Ply Surgical Face Mask Structure & Function



OUTER LAYER


Hydrophobic non-woven layer effectively blocks saliva

MIDDLE LAYER

Meltblown cloth filter non-oily particles in the air

INNER LAYER

Inner non-woven layer absorbs exhaled heat and leaves skin dry



Stallion Medical Inc.

Address: 3085 Beyer Blvd. Suite A-105, San Diego, CA. 92154

Phone (619) 428-6320

E-mail: sales@stallionmed.com

Follow us:



www.stallionmedical.com

We search the globe for the best materials