



Stallion Medical Inc.

ACP-1500

Auto Chart Projector



Specification

Projection distance: 1.5m ~ 7m

Projection magnification: 30 x (at 5m)

Projection size: 330mm(W) x 270mm(H) (at 5m)

Chart Switching speed: 0.01~0.02S

Mask Number: 1 open, 5 horizontal lines, 8 vertical lines, 21 single letters, 1 red/green

Mask Switching speed: One mask per 0.01~0.02s

Program: 2 sets programs, each program contains up to 30 steps

Lamp: 12V, 50W (halogen lamp) 12V, 10W(LED)


Power saving fuction: Will automatically standby in 10 min

Power source: AC 220V, 50Hz or 110V, 60Hz

Power consumption: 80W

Dimensions: 320mm(L) x 218mm(W) x 240mm(H)

Net weight: 7.2 Kg



Standard accesories

Remote control, Polarized metal screen, halogen lamp, polarized glasses, fuses(2), batteries(2)

Optional accesories

Floor stand, wall stand



Stallion Medical Inc.

Address: 3085 Beyer Blvd. Suite A-105, San Diego, CA. 92154

Phone (619) 428-6320 • Fax (619) 428-6323

E-mail: sales@stallionmed.com

Follow us:



www.stallion-medical.com