



Stallion Medical Inc.

Face Mask

KN95



Specification

Material: Non woven and meltblown grade 95

Size: 65 x 55 x 150mm

Layers: 4 Layers

FE (filter efficiency %): $\geq 95\%$

Inhaling Resistance [Pa]: Inhaling Resistance ≤ 350

Exhaling Resistance [Pa]: Exhaling Resistance ≤ 250


Executive Standard: GB2626-2006 (30 \pm 10)%

Storage: Temperature: (25 \pm 5) °C Humidity: (30 \pm 10)%

Validity: 2 Years

4 LAYERS OF PROTECTION

Filtration rate 95%



Stallion Medical Inc.

Address: 3085 Beyer Blvd. Suite A-105, San Diego, CA.

92154 Phone (619) 428-6320

E-mail: sales@stallionmed.com

Follow us:



www.stallionmedical.com

We search the globe for the best materials