



Stallion Medical Inc.

LYS-IB

Manual Lensmeter



Specification

Ranges of diopter measurement: 0~+25D

Minimum scale value: 0.125D at 0~±5D

0.25D at ±5D~±25D

Astigmatism axis of cylindrical lens 0~180°,Minimum scale value 5°

Prism diopter:

prism diopter of patten A:at 1Δinterval in the range of 0~5Δ

prism diopter of patten B:at 1Δinterval in the range of 0~20√

Prism base angle: 0~180°,Minimum scale value 5°

Regulation of ocular visibility: ±5D

Size of lens: φ16mm~ φ80mm

Overall dimensions of device: 275mm×130mm×455mm

Weight: 5.6kg

Lighting lamp: 220v15w 110v15w



Stallion Medical Inc.

Address: 3085 Beyer Blvd. Suite A-105, San Diego, CA.

92154 Phone (619) 428-6320

E-mail: sales@stallionmed.com

Follow us:



www.stallionmedical.com