



Stallion Medical Inc.

S350 Slit Lamp



Microscope

Microscope type: Galilean Type

Magnification Change: Revolving Drum(5 or 3 steps)**Total Magnification:** 10X, 16X

Total Magnification: 6X, 10X, 16X, 25X, 40X (S360 / S350) ; 10X, 16X, 25X (S360S / S350S)

Eyepieces: 12.5X

Angle between eyepieces: 13°

Pupillary Adjustment: 52mm - 78mm

Diopter Adjustment: ±6D

Field of View: Ø34.7mm, Ø22mm, Ø13.5mm, Ø8.5mm, Ø5.5mm(S360 / S350),

Ø22mm, Ø13.5mm, Ø8.5mm (S360S/ S350S)

Slit Illumination

Slit Width: 0~14mm continuous (slit becomes a circle at 14mm)

Slit Length: 1~14mm continuous

Aperture Diameters: ø14mm, ø10mm, ø5mm, ø3mm, ø2mm, ø1mm, ø0.2mm

Slit Angle: 0°~180°

Slit inclination: 5°, 10°, 15°, 20°

Filters: Heat-absorbing filter, ND filter, Red-Free filter, Cobalt Blue filter

Lamp: 6V/20W Halogen Lamp (S350 / S350S)

Luminance: ≥50KLX (S350 / S350S)

Instrument Base

Longitudinal Movement: 110mm

Lateral Movement: 110mm

Vertical Movement: 30mm

Fine Base Movement: 15mm

Chin Rest

Vertical Movement: 80mm

Fixation: LED

Power Supply

Input Voltage: 220V/110V ~ +/-10%

Input Frequency: 50Hz/60Hz

Power Consumption: 30VA (max)

Output Voltage: Lamp 6V, Fixation 3V

Packaging

Dimension: 770mm x 470mm x 570mm (L/W/H)

Gross weight: 24kg

Net weight: 16kg



Stallion Medical Inc.

Address: 3085 Beyer Blvd. Suite A-105, San Diego, CA. 92154

Phone (619) 428-6320 • Fax (619) 428-6323

E-mail: sales@stallionmed.com

Follow us:



www.stallion-medical.com