




Stallion Medical Inc.

# S360 / S360S

Slit Lamp



## Features

---

**Microscope Type:** Galilean Type

**Maginification Chance:** Revolving Drum (5 or 3 steps)

**Total Maginification:** 6X, 10X, 16X, 25X, 40X (S360); 10X, 16X, 25X (S360S)

**Eyepieces:** 12.5X

**Angle between Eyepieces:** 13°

**Pupillary Adjustment:** 52mm-78mm

**Diopter Adjustment:** +/- 6D

**Field of View:** ø34.7mm, ø22mm, ø13.5mm, ø8.5mm, ø5.5mm (S360), ø22mm, ø13.5mm, ø8.5 (S360S)

## Slit Illumination

---

**Slit Width:** 0~14mm

**Slit Length:** 1~14mm continuous

**Aperture Diameters:** ø14mm, ø10mm, ø5mm, ø3mm, ø2mm, ø1mm, ø0.2mm

**Slit Angle:** 0°~180°

**Slit Inclination:** 5°, 10°, 15°, 20°

**Filters:** Heat-absorbing filter, ND Filter, Red-free filter, Cobalt Blue Filter

**Lamp:** 3V LED Module (S360/S360S)

**Luminance:** ≥ 60KLX (S360/S360S)

## Instrument Base

---

**Longitudinal Movement:** 110mm

**Lateral Movement:** 110mm

**Vertical Movement:** 30mm

**Fine Base Movement:** 15mm

## Chin Rest

---

**Vertical Movement:** 80mm

**Fixation:** LED

## Power Supply

---

**Input Voltage:** 220V/220V~+/-10%

**Input Frequency:** 50Hz/60Hz

**Power Consumption:** 10VA (max) (S360/S360S)

**Output Voltage:** Lamp 6V, Fixation 3V

## Packaging

---

**Dimension:** 770mm x 470mm x 570mm (L/W/H)

**Gross weight:** 24kg

**Net weight:** 16kg



**Stallion Medical Inc.**

Address: 3085 Beyer Blvd. Suite A-105, San Diego, CA. 92154

Phone (619) 428-6320 • Fax (619) 428-6323

E-mail: sales@stallionmed.com

Síguenos en:



[www.stallion-medical.com](http://www.stallion-medical.com)